

---

# **paws Documentation**

***Release 0.4.0***

**Lenson A. Pellouchoud**

**Feb 17, 2017**



# CONTENTS

<b>1</b>	<b>Introduction</b>	<b>3</b>
<b>2</b>	<b>Quick Start</b>	<b>5</b>
<b>3</b>	<b>Installation</b>	<b>7</b>
3.1	Installing with pip . . . . .	7
3.2	Downloading Source . . . . .	7
3.3	Testing . . . . .	7
<b>4</b>	<b>API Documentation</b>	<b>9</b>
4.1	paws package . . . . .	9
<b>5</b>	<b>Indices and tables</b>	<b>41</b>
	<b>Python Module Index</b>	<b>43</b>



Contents:



## INTRODUCTION

The `paws` package aims to provide a fast and lean platform for building and executing workflows for data processing. It was originally developed to perform analysis of diffraction images for research purposes at SLAC/SSRL. At the core of `paws` is a workflow engine that uses a library of operations to crunch through data and expose select results while attempting to minimize resource consumption.

`paws` is currently written in Python, based on Qt via the PySide bindings. Internally, `paws` keeps track of data in Qt-based tree models, which can be controlled either directly (through the `paws` api) or through a gui (employing the Qt model-view framework).

`paws` also provides an interface to `xi-cam`, a synchrotron x-ray diffraction data analysis package written by the CAMERA Institute and Pandolfi, et al at the Lawrence Berkeley National Lab.

Some the core goals of `paws`:

- Eliminate redundant development efforts
- Streamline and standardize routine data analysis
- Simplify data storage and provide large-scale analysis
- Perform data analysis in real time for results-driven feedback

The `paws` developers would love to hear from you if you have wisdom, thoughts, haikus, bugs, artwork, or suggestions. Limericks are also welcome. Get in touch with us at `paws-developers@slac.stanford.edu`.





## QUICK START

Minimal and usually-effective installation instructions.

Here is a reference to the *[brief introduction](#)*.

*This chapter* is for setting up `paws` quickly in an environment that is prepared to install Python packages with `pip`.



## INSTALLATION

Here are instructions for installing `paws` from PyPI, or downloading and testing the `paws` source code.

### Installing with pip

Instructions will go here for installing `paws` using the Python package installer `pip` (currently not implemented).

### Downloading Source

The source code for `paws` is hosted on `github`. Clone the repository from `https://github.com/slaclab/paws.git`. You should then be able to run `paws` by invoking `python main.py` from the root directory.

### Testing

`paws` comes with a tests that can be used to ensure the platform runs as expected. After *downloading the source*, invoke `python -m unittest discover` from the root directory.



## API DOCUMENTATION

This is the complete auto-generated documentation of the `paws` package, made with sphinx-apidoc.

### `paws` package

#### Subpackages

##### `paws.core` package

##### Subpackages

##### `paws.core.operations` package

##### Subpackages

##### `paws.core.operations.COMM` package

##### Submodules

##### `paws.core.operations.COMM.TestTCP` module

Operation for testing a TCP client

```
class paws.core.operations.COMM.TestTCP.TestTCP
```

Bases: `paws.core.operations.operation.Operation`

Given a filename and a TCP client, this operation sends some stuff to the TCP client.

```
run ()
```

##### Module contents

##### `paws.core.operations.DISPLAY` package

##### Submodules

### paws.core.operations.DISPLAY.simple\_plots module

```
class paws.core.operations.DISPLAY.simple_plots.MPLFigFromXYData
    Bases: paws.core.operations.operation.Operation

    run ()

class paws.core.operations.DISPLAY.simple_plots.SimplePlot
    Bases: paws.core.operations.operation.Operation

    Plot a 1d array against another 1d array.

    run ()

    simple_plot (x, y)
```

### paws.core.operations.DISPLAY.zip module

```
class paws.core.operations.DISPLAY.zip.LogLogZip
    Bases: paws.core.operations.operation.Operation

    Takes the logarithm of two 1d ndarrays, then zips them together.

    Logarithm is taken in base ten.

    Any elements with non-positive values are removed, so this operation may not be appropriate for computational
    purposes.

    run ()

class paws.core.operations.DISPLAY.zip.Zip
    Bases: paws.core.operations.operation.Operation

    Zips two 1d ndarrays together.

    run ()

paws.core.operations.DISPLAY.zip.zip (x, y)
    Zips input 1d vectors together for display.

    Should be pre-checked.

paws.core.operations.DISPLAY.zip.zip_check (x, y)
    Checks that inputs are 1d vectors of the same size and shape.
```

## Module contents

### paws.core.operations.EXECUTION package

#### Submodules

### paws.core.operations.EXECUTION.batch\_execution module

```
class paws.core.operations.EXECUTION.batch_execution.BatchFromFiles
    Bases: paws.core.operations.operation.Batch

    Provides a sequence of inputs to be used in repeated execution of a workflow. Collects the outputs produced for
    each of the inputs.
```

**batch\_ops()**  
Provide a list of uri's of ops to be included in batch execution

**input\_list()**

**input\_routes()**  
Provide the input route in a list- must return list.

**output\_list()**

**run()**  
Build a list of [uri:value] dicts to be used in the workflow.

**saved\_items()**  
List uris to be saved/stored after execution

### **paws.core.operations.EXECUTION.realtime\_execution module**

**class** `paws.core.operations.EXECUTION.realtime_execution.RealtimeFromFiles`  
Bases: `paws.core.operations.operation.Realtime`

Provides inputs to be used in repeated execution of a workflow from files with names matching a regex, as they arrive in a specified directory. Collects the outputs produced for each of the inputs.

**static delay()**  
Amount of time to wait between execution attempts, in milliseconds

**input\_iter()**

**input\_routes()**  
Use the Realtime.input\_locators to list uri's of all input routes- must return list.

**output\_list()**

**realtime\_ops()**  
Use the Realtime.input\_locator to list uri's of ops to be saved/stored after execution

**run()**  
This should create an iterator whose next() gives a {uri:value} dict built from the latest-arrived file

**saved\_items()**  
Use the Realtime.input\_locator to list uri's of ops to be included in realtime execution

### **Module contents**

#### **paws.core.operations.INPUT package**

#### **Subpackages**

#### **paws.core.operations.INPUT.BL1-5 package**

#### **Submodules**

#### **paws.core.operations.INPUT.BL1-5.SSRL\_1\_5\_readers module**

### **Module contents**

## paws.core.operations.INPUT.WXDIFF package

### Submodules

#### paws.core.operations.INPUT.WXDIFF.read\_wxdiff\_files module

**class** `paws.core.operations.INPUT.WXDIFF.read_wxdiff_files.WXDCalibToDict`  
Bases: `paws.core.operations.operation.Operation`  
Input is the path to a WXDiff .calib file output is a dict containing calib parameters  
**run** ()

### Module contents

### Submodules

#### paws.core.operations.INPUT.fabio\_input module

**class** `paws.core.operations.INPUT.fabio_input.LoadFabIO`  
Bases: `paws.core.operations.operation.Operation`  
Takes a filesystem path and calls fabIO to load it.  
**run** ()  
Call on fabIO to extract image data

#### paws.core.operations.INPUT.pillow\_tif\_input module

**class** `paws.core.operations.INPUT.pillow_tif_input.LoadTif`  
Bases: `paws.core.operations.operation.Operation`  
Takes a filesystem path that points to a .tif, outputs image data and metadata from the file.  
**istiff** ()  
**run** ()  
Call on image rendering libs to extract image data \*.tiff or \*.tif images: use PIL  
**tifftest** = `<_sre.SRE_Pattern object>`

**class** `paws.core.operations.INPUT.pillow_tif_input.LoadTifToGrayscale`  
Bases: `paws.core.operations.operation.Operation`  
Takes a filesystem path that points to a .tif, outputs image data and metadata from the file.  
**istiff** ()  
**run** ()  
Call on image rendering libs to extract image data \*.tiff or \*.tif images: use PIL  
**tifftest** = `<_sre.SRE_Pattern object>`



## paws.core.operations.INPUT.read\_csv module

**class** `paws.core.operations.INPUT.read_csv.ReadCSV_q_I_dI`

Bases: `paws.core.operations.operation.Operation`

Read q, I, and (if available) dI from a csv-formatted file.

If the csv has no third column, returns None for dI.

**run** ()

`paws.core.operations.INPUT.read_csv.read_csv_q_I_maybe_dI` (*nameloc*)

## Module contents

### paws.core.operations.OUTPUT package

#### Subpackages

#### paws.core.operations.OUTPUT.PIF package

#### Submodules

### paws.core.operations.OUTPUT.PIF.citration\_shipment module

**class** `paws.core.operations.OUTPUT.PIF.citration_shipment.CheckDataSet`

Bases: `paws.core.operations.operation.Operation`

Take a Citration client as input and use it to query a data set. Output some indication of whether or not the query was successful.

**run** ()

**class** `paws.core.operations.OUTPUT.PIF.citration_shipment.CreateDataSet`

Bases: `paws.core.operations.operation.Operation`

Take a Citration client as input and use it to create a data set. Output the index of the created data set.

**run** ()

**class** `paws.core.operations.OUTPUT.PIF.citration_shipment.ShipJSON`

Bases: `paws.core.operations.operation.Operation`

Take a .json file containing a pif, and ship the pif to a given Citration data set.

**run** ()

**class** `paws.core.operations.OUTPUT.PIF.citration_shipment.ShipToDataSet`

Bases: `paws.core.operations.operation.Operation`

Take a `pycif.obj.System` object and ship it to a given Citration data set.

**run** ()

## Module contents

### Submodules

#### paws.core.operations.OUTPUT.write\_csv module

**class** `paws.core.operations.OUTPUT.write_csv.WriteArrayCSV`

Bases: `paws.core.operations.operation.Operation`

Write a 2d array to a csv file

**run** ()

**class** `paws.core.operations.OUTPUT.write_csv.WriteCSV_q_I_dI`

Bases: `paws.core.operations.operation.Operation`

Write q, I, and (if available) dI to a csv-formatted file.

**run** ()

`paws.core.operations.OUTPUT.write_csv.replace_extension (old_name, new_extension)`

Return a file name that is identical except for extension.

#### Parameters

- **old\_name** – string path or file name
- **new\_extension** – string extension, e.g. ".txt"

#### Returns

Accepts extensions with or without an initial ".".

`paws.core.operations.OUTPUT.write_csv.write_csv_q_I_maybe_dI (q, I, dI, nameloc)`

## Module contents

### paws.core.operations.PACKAGING package

#### Subpackages

#### paws.core.operations.PACKAGING.PIF package

### Submodules

#### paws.core.operations.PACKAGING.PIF.pif\_saxs module

**class** `paws.core.operations.PACKAGING.PIF.pif_saxs.PifNPSynth`

Bases: `paws.core.operations.operation.Operation`

Package results from nanoparticle solution synthesis into a pypif.obj.ChemicalSystem object.

**run** ()

**saxs\_to\_pif\_properties** (*q, I<sub>q</sub>, T<sub>C</sub>*)

## paws.core.operations.PACKAGING.PIF.pif\_ternary module

```
class paws.core.operations.PACKAGING.PIF.pif_ternary.PifTernary
    Bases: paws.core.operations.operation.Operation

    Package results from ternary wafer synthesis into a pypif.obj.ChemicalSystem object.

    intensity_features_to_pif_properties (I_feats)

    pks_to_pif_properties (q_pk, I_pk)

    run ()

    texture_to_pif_properties (q, tex_q)

    xrd_to_pif_properties (q, I_q)
```

## Module contents

### Submodules

## paws.core.operations.PACKAGING.data\_packing module

```
class paws.core.operations.PACKAGING.data_packing.TimeTempFromHeader
    Bases: paws.core.operations.operation.Operation

    Get time and temperature from a detector output header file. Return string time, float time (utc in seconds), and float temperature. Time is assumed to be in the format Day Mon dd hh:mm:ss yyyy.

    run ()

class paws.core.operations.PACKAGING.data_packing.WindowZip
    Bases: paws.core.operations.operation.Operation

    From input iterables of x and y, produce an n-by-2 array where x is bounded by the specified limits

    run ()

class paws.core.operations.PACKAGING.data_packing.Window_q_I_2
    Bases: paws.core.operations.operation.Operation

    From input iterables of q_list_in and I_list_in, produce two 1d vectors where q is greater than q_min and less than q_min

    run ()

class paws.core.operations.PACKAGING.data_packing.XYDataFromBatch
    Bases: paws.core.operations.operation.Operation

    Given a batch output and appropriate keys, use the uris to harvest x and y data from the batch.

    run ()
```

## Module contents

## paws.core.operations.PROCESSING package

### Subpackages

## paws.core.operations.PROCESSING.SAXS package

### Submodules

#### paws.core.operations.PROCESSING.SAXS.background\_subtraction module

**class** `paws.core.operations.PROCESSING.SAXS.background_subtraction.BgSubtractByTemperature`  
Bases: `paws.core.operations.operation.Operation`

Find a background spectrum from a batch of background spectra, where the temperature of the background spectrum is as close as possible to the (input) temperature of the measured spectrum. Then subtract that background spectrum from the input spectrum. The measured and background spectra are expected to have the same domain.

**run** ()

**class** `paws.core.operations.PROCESSING.SAXS.background_subtraction.SubtractMaximumBackground`  
Bases: `paws.core.operations.operation.Operation`

Subtract a background from a foreground, with scaling to prevent over-subtraction.

Has optional arguments for error vectors. If no error estimate is available, set these to *None*. If either error vector is *None*, the output error vector will also be *None*.

Can also, optionally, window to a region in *q*.

**run** ()

`paws.core.operations.PROCESSING.SAXS.background_subtraction.subtract_maximum_background` (*q\_f*,  
*q\_b*,  
*I\_f*,  
*I\_b*,  
*dI\_f*,  
*dI\_b*,  
*q1*=  
*q2*=

#### paws.core.operations.PROCESSING.SAXS.saxs module

**class** `paws.core.operations.PROCESSING.SAXS.saxs.FetchReferences`  
Bases: `paws.core.operations.operation.Operation`

Fetch previously generated and stored metrics used for guessing diffraction pattern properties.

This function is for examination, and is not necessary for computation. The file in question is auto-fetched by operations that rely on it.

**run** ()

**class** `paws.core.operations.PROCESSING.SAXS.saxs.GenerateReferences`  
Bases: `paws.core.operations.operation.Operation`

Generate metrics used for guessing diffraction pattern properties.

DOES NOT STORE. IS NOT AUTOMATICALLY AVAILABLE TO OTHER OPERATIONS.

Use `OverwriteReferences` to store the resulting dictionary, if desired.

**run** ()

**class** `paws.core.operations.PROCESSING.SAXS.saxs.GenerateSphericalDiffraction`  
 Bases: `paws.core.operations.operation.Operation`

Generate a SAXS diffraction pattern for spherical nanoparticles.

Uses `r0` in real units and, if asked to generate `q`, gives `q` in real units.

**run** ()

**class** `paws.core.operations.PROCESSING.SAXS.saxs.GuessProperties`  
 Bases: `paws.core.operations.operation.Operation`

Guess the polydispersity, mean size, and amplitude of spherical diffraction pattern.

Assumes the data have already been background subtracted, smoothed, and otherwise cleaned.

**run** ()

**class** `paws.core.operations.PROCESSING.SAXS.saxs.OptimizeSphericalDiffractionFit`  
 Bases: `paws.core.operations.operation.Operation`

From an initial guess, optimize `r0`, `I0`, and `fractional_variation`.

**run** ()

**class** `paws.core.operations.PROCESSING.SAXS.saxs.OverwriteReferences`  
 Bases: `paws.core.operations.operation.Operation`

Store previously generated metrics used for guessing diffraction pattern properties.

USE WITH CAUTION as it will OVERWRITE any existing reference file.

**run** ()

`paws.core.operations.PROCESSING.SAXS.saxs.arbitrary_order_solution` (*order*, *x*, *y*,  
*dy=None*)

Solves for a polynomial “fit” of arbitrary order.

`paws.core.operations.PROCESSING.SAXS.saxs.blur` (*x*, *factor*)

`paws.core.operations.PROCESSING.SAXS.saxs.chi_squared` (*y1*, *y2*)

`paws.core.operations.PROCESSING.SAXS.saxs.choose_dips_and_shoulders` (*q*, *I*,  
*dI=None*)

Find the location of dips (low points) and shoulders (high points).

`paws.core.operations.PROCESSING.SAXS.saxs.clean_extrema` (*dips*, *shoulders*, *y*)

`paws.core.operations.PROCESSING.SAXS.saxs.downwards_weaksauce_identifier` (*four\_indices*,  
*y*)

`paws.core.operations.PROCESSING.SAXS.saxs.dummy` (*array1d*)

Turn 1d array into dummy-index vector for 2d matrix computation.

Sum over the dummy index by taking *object.sum(axis=0)*.

`paws.core.operations.PROCESSING.SAXS.saxs.dump_references` (*references*)

`paws.core.operations.PROCESSING.SAXS.saxs.factorStretched` (*x*, *factor*)

`paws.core.operations.PROCESSING.SAXS.saxs.first_dip` (*q*, *I*, *dips*, *dI=None*)

`paws.core.operations.PROCESSING.SAXS.saxs.fullFunction` (*x*)

`paws.core.operations.PROCESSING.SAXS.saxs.gauss` (*x*, *x0*, *sigma*)

`paws.core.operations.PROCESSING.SAXS.saxs.gauss_guess` (*signalMagnitude*, *signalCurvature*)

Guesses a gaussian intensity and width from signal magnitude and curvature.

**Parameters**

- **signalMagnitude** – number-like with units of magnitude
- **signalCurvature** – number-like with units of magnitude per distance squared

**Return intensity, sigma**

The solution given is not fitted; it is a first estimate to be used in fitting.

`paws.core.operations.PROCESSING.SAXS.saxs.gen_q_vector(qmin, qmax, qstep)`  
An inclusive-range version of `np.arange`.

`paws.core.operations.PROCESSING.SAXS.saxs.generateRhoFactor(factor)`  
Generate a gaussian distribution of number densities (`rho`).

**Parameters** `factor` – float

**Returns**

`factor` should be 0.1 for a sigma 10% of size

`paws.core.operations.PROCESSING.SAXS.saxs.generate_references(x, factorVals)`

`paws.core.operations.PROCESSING.SAXS.saxs.generate_spherical_diffraction(q,  
i0,  
r0,  
poly)`

`paws.core.operations.PROCESSING.SAXS.saxs.guess_nearest_point_on_nonmonotonic_trace_normal`

Finds the nearest point to location *loclist* on a trace *tracelist*.

Uses normalized distances, i.e., the significance of each dimension is made equivalent.

*loclist* and *tracelist* are lists of the same length and type. Elements of *loclist* are floats; elements of *tracelist* are arrays of the same shape as *coordinate*. *coordinate* is the independent variable along which *tracelist* is varying.

*tracelist* must have decent sampling to start with or the algorithm will likely fail.

`paws.core.operations.PROCESSING.SAXS.saxs.guess_noise_floor(q, I, r0)`

`paws.core.operations.PROCESSING.SAXS.saxs.guess_polydispersity(q, I, dI=None)`

`paws.core.operations.PROCESSING.SAXS.saxs.guess_size(fractional_variation,  
first_dip_q)`

`paws.core.operations.PROCESSING.SAXS.saxs.horizontal(array1d)`

Turn 1d array into 2d horizontal vector.

`paws.core.operations.PROCESSING.SAXS.saxs.load_references()`

`paws.core.operations.PROCESSING.SAXS.saxs.local_maxima_detector(y)`

Finds local maxima in ordered data *y*.

**Parameters** `y` – 1d numpy float array

**Return maxima** 1d numpy bool array

*maxima* is *True* at a local maximum, *False* otherwise.

This function makes no attempt to reject spurious maxima of various sorts. That task is left to other functions.

`paws.core.operations.PROCESSING.SAXS.saxs.local_minima_detector(y)`

Finds local minima in ordered data *y*.

**Parameters** *y* – 1d numpy float array

**Return minima** 1d numpy bool array

*minima* is *True* at a local minimum, *False* otherwise.

This function makes no attempt to reject spurious minima of various sorts. That task is left to other functions.

`paws.core.operations.PROCESSING.SAXS.saxs.make_poly_matrices(x, y, error, order)`

Make the matrices necessary to solve a polynomial fit of order *order*.

**Parameters**

- *x* – 1d array representing independent variable
- *y* – 1d array representing dependent variable
- *error* – 1d array representing uncertainty in *y*
- *order* – integer order of polynomial fit

**Return matrix, vector** MC=V, where M is *matrix*, V is *vector*, and C is the polynomial coefficients to be solved for. *matrix* is an array of shape (order+1, order+1). *vector* is an array of shape (order+1, 1).

`paws.core.operations.PROCESSING.SAXS.saxs.polydispersity_metric_heightAtZero(qFirstDip, q, I, dI=None)`

`paws.core.operations.PROCESSING.SAXS.saxs.polydispersity_metric_sigmaScaledFirstDip(q, I, dips, shoulders, qFirstDip, heightFirstDip, dI=None)`

`paws.core.operations.PROCESSING.SAXS.saxs.polynomial_value(coefficients, x)`

Finds the value of a polynomial at a location.

`paws.core.operations.PROCESSING.SAXS.saxs.power_law_solution(x, y, dy=None)`

Solves for a power law by solving for a linear fit in log-log space.

`paws.core.operations.PROCESSING.SAXS.saxs.prep_for_pickle(factorVals, xFirstDip, sigmaScaledFirstDip, heightFirstDip, heightAtZero)`

`paws.core.operations.PROCESSING.SAXS.saxs.quadratic_extremum(coefficients)`

Finds the location in independent coordinate of the extremum of a quadratic equation.

`paws.core.operations.PROCESSING.SAXS.saxs.refine_guess(q, I, I0, r0, frac, q1, I1)`

`paws.core.operations.PROCESSING.SAXS.saxs.take_polydispersity_metrics(x, y, dy=None)`

```
paws.core.operations.PROCESSING.SAXS.saxs.upwards_weaksauce_identifier (four_indices,
                                                                    y)

paws.core.operations.PROCESSING.SAXS.saxs.vertical (array1d)
    Turn 1d array into 2d vertical vector.
```

## paws.core.operations.PROCESSING.SAXS.saxs\_features module

```
class paws.core.operations.PROCESSING.SAXS.saxs_features.NanoparticleFeatures
    Bases: paws.core.operations.operation.Operation

    Use a saxs spectrum (I(q) vs. q) to extract key features of a nanoparticle colloid: average size, standard deviation
    of sizes, and nanoparticle density. This module assumes spherical nanoparticles.

    static monodisperse_saxs (q, r)
        Generate spherical nanoparticle saxs spectrum for monodisperse particle radius r at scattering vectors q.

    static polydisperse_saxs (q, r0, p)
        Generate spherical nanoparticle saxs spectrum for mean particle radius r0, size polydispersity (sigma(r)/r0)
        p, and scattering vectors q.

    run ()
```

## Module contents

### Submodules

## paws.core.operations.PROCESSING.HiTp\_feature\_extraction module

Created on Mon Jun 06 18:02:32 2016

author: Fang Ren, Apurva Mehta For details, refer to the recent paper submitted to ACS Combinatorial Science Fang Ren implemented all the methods on scripting level, and originally contributed to slacx

Contributor: fangren (please add if there is more)

```
class paws.core.operations.PROCESSING.HiTp_feature_extraction.FindPeaksByWindow
    Bases: paws.core.operations.operation.Operation

    Walk a 1d array and find its local maxima. A maximum is found if it is the highest point within window size of
    itself. An optional threshold for the peak intensity relative to the window-average can be used to filter out peaks
    due to noise.

    run ()

class paws.core.operations.PROCESSING.HiTp_feature_extraction.IntensityFeatures
    Bases: paws.core.operations.operation.Operation

    Extract the maximum intensity, average intensity, and a ratio of the two from data

    run ()

class paws.core.operations.PROCESSING.HiTp_feature_extraction.PeakFeatures
    Bases: paws.core.operations.operation.Operation

    Extract the locations and intensities of peaks from a 1D spectrum

    get_extrema (x, y, delta)
        Given vectors x and y, return an n-by-2 array of x,y pairs for the minima and maxima of y(x)
```



```
run ()
```

```
class paws.core.operations.PROCESSING.HitP_feature_extraction.TextureFeatures
    Bases: paws.core.operations.operation.Operation
```

Analyze the texture

```
run ()
```

### paws.core.operations.PROCESSING.VoigtPeakFit module

```
class paws.core.operations.PROCESSING.VoigtPeakFit.VoigtPeakFit
    Bases: paws.core.operations.operation.Operation
```

Fit a set of x and y values to a Voigt distribution. Solves the distribution over the space of hwhm (half width at half max) of the gaussian and lorentzian distributions and distribution center. Takes as input a guess for the distribution center and hwhm. Range of fit is determined by weighting the objective by a Hann window centered at the distribution center, with a window width of the distribution's estimated full width at half max.

```
static gaussian (x, hwhm_g)
    gaussian distribution at points x, center 0, hwhm hwhm_g
```

```
static hann_voigt_fit (x, y, xc, hwhm_g, hwhm_l, scl)
```

```
static lorentzian (x, hwhm_l)
    lorentzian distribution at points x, center 0, hwhm hwhm_l
```

```
run ()
```

```
static solve_voigt (x, y, xc, hwhm_g, hwhm_l, scl)
    iteratively minimize an objective to fit x, y curve to a voigt profile
```

```
static voigt (x, hwhm_g, hwhm_l)
    voigt distribution resulting from convolution of a gaussian with hwhm hwhm_g and a lorentzian with
    hwhm hwhm_l
```

### paws.core.operations.PROCESSING.arithmetic module

```
class paws.core.operations.PROCESSING.arithmetic.Add
    Bases: paws.core.operations.operation.Operation
```

Add two objects.

```
run ()
```

```
class paws.core.operations.PROCESSING.arithmetic.Divide
    Bases: paws.core.operations.operation.Operation
```

Divide two objects.

```
run ()
```

```
class paws.core.operations.PROCESSING.arithmetic.Exponentiate
    Bases: paws.core.operations.operation.Operation
```

Exponentiate an object by another object.

```
run ()
```

**class** `paws.core.operations.PROCESSING.arithmetic.Logarithm`

Bases: `paws.core.operations.operation.Operation`

Take the logarithm of an object by some base.

**run** ()

**class** `paws.core.operations.PROCESSING.arithmetic.Multiply`

Bases: `paws.core.operations.operation.Operation`

Multiply two objects.

**run** ()

**class** `paws.core.operations.PROCESSING.arithmetic.Subtract`

Bases: `paws.core.operations.operation.Operation`

Subtract one object from another.

**run** ()

## **paws.core.operations.PROCESSING.calib\_and\_reduce module**

Operations for remeshing and and reducing an image

author: Fang Ren, Apurva Mehta, Ron Pandolfi For details, please refer to the recent paper submitted to ACS Combinatorial Science F. Ren implemented the remeshing on scripting level, contributed by R. Pandolfi from LBL. F. Ren originally contributed it to slax

Please add yours too, Amanda, Lenson.

contributors: fangren, apf, lensonp

**class** `paws.core.operations.PROCESSING.calib_and_reduce.CalByWXDDict`

Bases: `paws.core.operations.operation.Operation`

Input image data (ndarray) and a dict of calibration parameters from a WxDiff .calib file Return q, chi, I(q,chi)

**run** ()

**class** `paws.core.operations.PROCESSING.calib_and_reduce.ReduceByWXDDict`

Bases: `paws.core.operations.operation.Operation`

Input image data (ndarray) and a dict of calibration parameters from a WxDiff .calib file Return q, I(q)

**run** ()

**class** `paws.core.operations.PROCESSING.calib_and_reduce.ReduceByWXDDict_mask_error`

Bases: `paws.core.operations.operation.Operation`

Input image data (ndarray) and a dict of calibration parameters from a WxDiff .calib file Return q, I(q)

**run** ()

## **paws.core.operations.PROCESSING.mirror module**

**class** `paws.core.operations.PROCESSING.mirror.ArrayMirrorHorizontal`

Bases: `paws.core.operations.operation.Operation`

Mirror an array across a horizontal plane, i.e., exchange indices along axis 0.

**run** ()

**class** `paws.core.operations.PROCESSING.mirror.ArrayMirrorVertical`  
 Bases: `paws.core.operations.operation.Operation`

Mirror an array across a vertical plane, i.e., exchange indices along axis 1.

**run** ()

**class** `paws.core.operations.PROCESSING.mirror.MirrorHorizontally`  
 Bases: `paws.core.operations.operation.Operation`

Mirror an image, exchanging left and right.

I.e., mirror an ndarray along axis = 0.

**class** `paws.core.operations.PROCESSING.mirror.MirrorVertically`  
 Bases: `paws.core.operations.operation.Operation`

Mirror an image, exchanging top and bottom.

I.e., mirror an ndarray along axis = 1.

**run** ()

### **paws.core.operations.PROCESSING.rotation module**

**class** `paws.core.operations.PROCESSING.rotation.Rotation`  
 Bases: `paws.core.operations.operation.Operation`

Rotate an array by 90, 180, or 270 degrees.

**run** ()

Rotate self.inputs['image\_data'] and save as self.outputs['image\_data']

### **paws.core.operations.PROCESSING.smoothing module**

**class** `paws.core.operations.PROCESSING.smoothing.RectangularSmooth`  
 Bases: `paws.core.operations.operation.Operation`

Applies rectangular (moving average) smoothing filter to 1d data.

User-specified error estimate used to weight points. Set *dy* to *None* if unavailable.

**run** ()

**class** `paws.core.operations.PROCESSING.smoothing.SavitzkyGolay`  
 Bases: `paws.core.operations.operation.Operation`

Applies a Savitzky-Golay (polynomial fit approximation) filter to 1d data.

Uses error bars on intensity if available. Set *dy* to *None* otherwise.

**run** ()

**class** `paws.core.operations.PROCESSING.smoothing.TriangularSmooth`  
 Bases: `paws.core.operations.operation.Operation`

Applies triangular-weighted (moving average) smoothing filter to 1d data.

User-specified error estimate used to weight points. Set *dy* to *None* if unavailable.

**run** ()

`paws.core.operations.PROCESSING.smoothing.choose_m(order, base)`

Choose a number of points to use for SG fit.

#### Parameters

- **order** – integer order of polynomial fit
- **base** – an edge condition specification parameter

**Return m** integer number of data points to use

Helper function for *savitzky\_golay*.

`paws.core.operations.PROCESSING.smoothing.choose_start_and_end(m, base, ii, size)`

`paws.core.operations.PROCESSING.smoothing.dummy(array1d)`

Turn 1d array into dummy-index vector for 2d matrix computation.

Sum over the dummy index by taking *object.sum(axis=0)*.

`paws.core.operations.PROCESSING.smoothing.horizontal(array1d)`

Turn 1d array into 2d horizontal vector.

`paws.core.operations.PROCESSING.smoothing.make_poly_matrices(x, y, error, order)`

Make the matrices necessary to solve a polynomial fit of order *order*.

#### Parameters

- **x** – 1d array representing independent variable
- **y** – 1d array representing dependent variable
- **error** – 1d array representing uncertainty in y
- **order** – integer order of polynomial fit

**Return matrix, vector** MC=V, where M is *matrix*, V is *vector*, and C is the polynomial coefficients to be solved for. *matrix* is an array of shape (order+1, order+1). *vector* is an array of shape (order+1, 1).

`paws.core.operations.PROCESSING.smoothing.moving_average(data, m, shape, error=None)`

`paws.core.operations.PROCESSING.smoothing.no_specified_error_weights(data)`

`paws.core.operations.PROCESSING.smoothing.polynomial(x, coefficients)`

Evaluate a polynomial with given coefficients at location x.

#### Parameters

- **x** – numeric value of coordinate
- **coefficients** – 1d ndarray of one or more elements; zeroth element is zero-order coefficient, 1st is 1st-order coefficient, etc.

**Return value** value of the specified polynomial at x

In *coefficients*, zeroth element is zero-order coefficient (i.e. constant offset), 1st is 1st-order coefficient (i.e. linear slope), etc.

`paws.core.operations.PROCESSING.smoothing.savitzky_golay(x, y, order, base, dy=None)`

`paws.core.operations.PROCESSING.smoothing.specified_error_weights(error)`

`paws.core.operations.PROCESSING.smoothing.square_weighting(n)`

`paws.core.operations.PROCESSING.smoothing.triangular_weighting(n)`

`paws.core.operations.PROCESSING.smoothing.vertical` (*array1d*)  
Turn 1d array into 2d vertical vector.

## Module contents

### `paws.core.operations.TESTS` package

#### Submodules

#### `paws.core.operations.TESTS.identity` module

**class** `paws.core.operations.TESTS.identity.Identity`  
Bases: `paws.core.operations.operation.Operation`  
An Operation testing class, loads its input into its output  
**run** ()

## Module contents

#### Submodules

#### `paws.core.operations.op_manager` module

**class** `paws.core.operations.op_manager.OpManager` (*\*\*kwargs*)  
Bases: `paws.core.treemodel.TreeSelectionModel`  
Tree structure for categorized storage and retrieval of Operations.  
**add\_cat** (*new\_cat, parent*)  
Add a category to the tree under parent if not already there. Return its index.  
**add\_op** (*op, parent*)  
add op to the tree under QModelIndex parent  
**data** (*item\_indx, data\_role*)  
**get\_op** (*indx*)  
**get\_op\_byname** (*op\_name*)  
**headerData** (*section, orientation, data\_role*)  
**idx\_of\_cat** (*catname, parent*)  
If cat exists under parent, return its index, else return an invalid QModelIndex  
**list\_op\_names** ()  
**load\_cats** (*cat\_list*)  
**load\_ops** (*cat\_op\_list*)  
Load OpManager tree from input cat\_op\_list. Format of cat\_op\_list is [(category1,op1),(category2,op2),...]. i.e. each operation in cat\_op\_list is specified by a tuple, where the first element is a category, and the second element is the Operation itself. load\_cats() MUST be called before load\_ops() and MUST ensure that all cats in cat\_op\_list exist in the tree.  
**remove\_op** (*removal\_indx*)

```
save_config()

staticMetaObject = <PySide.QtCore.QMetaObject object>
```

### paws.core.operations.operation module

```
class paws.core.operations.operation.Batch(input_names, output_names)
    Bases: paws.core.operations.operation.Operation

    batch_ops()
        Return a list of operation uris to be included in the Batch execution stack.

    input_list()
        Produce a list of OrderedDicts representing each set of inputs for the Batch to run. Each OrderedDict
        should be populated with [input_uri:input_value] pairs.

    input_routes()
        Produce a list of the input routes used by the Batch, in the same order as each of the OrderedDicts provided
        by Batch.input_list()

    output_list()
        Produce a list of OrderedDicts representing the outputs for each batch input. Each OrderedDict should be
        populated with [input_uri:input_value] pairs.

    saved_items()
        Return a list of items to be saved after each execution.

class paws.core.operations.operation.Operation(input_names, output_names)
    Bases: object

    description()
        self.description() returns a string documenting the input and output structure and usage instructions for the
        Operation

    doc_as_string()

    input_description()

    load_defaults()

    output_description()

    run()
        Operation.run() should use all of the items in Operation.inputs and set values for all of the items in Oper-
        ation.outputs.

class paws.core.operations.operation.Realtime(input_names, output_names)
    Bases: paws.core.operations.operation.Operation

    delay()
        Return the number of MILLIseconds to pause between iterations. Overload this method to change the
        pause time- default is 1 second.

    input_iter()
        Produce an iterator over OrderedDicts representing each set of inputs to run. Each dict should be populated
        with [input_uri:input_value] pairs. When there is no new set of inputs to run, should return None.

    input_routes()
        Produce a list of [input_uri] routes in the same order as the OrderedDicts produced by Realtime.input_iter()
```

**output\_list()**

Produce a list of OrderedDicts representing the outputs for each realtime input. Each OrderedDict should be populated with [input\_uri:input\_value] pairs.

**realtime\_ops()**

Return a list of operation uris to be included in the Realtime execution stack.

**saved\_items()**

Return a list of item uris to be saved after each execution.

**paws.core.operations.optools module**

Operations config and processing routines

**class** `paws.core.operations.optools.InputLocator` (*src=0, tp=0, val=None*)

Bases: `object`

Objects of this class are used as containers for inputs to an Operation, and should by design contain the information needed to find the relevant input data. After the data is loaded, it should be stored in `InputLocator.data`.

`paws.core.operations.optools.cast_type_val` (*tp, val*)

Perform type casting for operation inputs. This should be called only for source = text\_input.

`paws.core.operations.optools.dict_contains_uri` (*uri, d*)

`paws.core.operations.optools.get_uri_from_dict` (*uri, d*)

`paws.core.operations.optools.op_dict` (*op*)

`paws.core.operations.optools.op_inputs_dict` (*op*)

`paws.core.operations.optools.parameter_doc` (*name, value, doc*)

`paws.core.operations.optools.plugin_dict` (*pgin*)

`paws.core.operations.optools.print_stack` (*stk*)

`paws.core.operations.optools.stack_contains` (*itm, stk*)

`paws.core.operations.optools.stack_size` (*stk*)

**Module contents**

`paws.core.operations.load_cfg` (*cfg\_file*)

`paws.core.operations.load_ops_from_module` (*mod, cat\_root*)

`paws.core.operations.load_ops_from_path` (*path\_, pkg, cat\_root='MISC'*)

`paws.core.operations.save_cfg` (*cfg\_data, cfg\_file*)

**paws.core.plugins package****Submodules****paws.core.plugins.CitrinationPlugin module**

**class** `paws.core.plugins.CitrinationPlugin.CitrinationPlugin`

Bases: `paws.core.plugins.plugin.PawsPlugin`

Wrapper contains a Citrination client and implements the PawsPlugin abc interface.

```
content ()
description ()
ship_dataset (pifs)
start ()
stop ()
```

### paws.core.plugins.SpecClientPlugin module

```
class paws.core.plugins.SpecClientPlugin.SpecClientPlugin
    Bases: paws.core.plugins.plugin.PawsPlugin
    content ()
    description ()
    receiveLine ()
    sendCmd (cmd)
    sendLine (line)
    send_commands (cmd_list)
    send_text (txt)
    start ()
    stop ()
```

### paws.core.plugins.TCPClientPlugin module

```
class paws.core.plugins.TCPClientPlugin.TCPClientFactory (protocol)
    Bases: twisted.internet.protocol.ClientFactory
    buildProtocol (addr)
    clientConnectionFailed (connector, reason)
        Clients call this when they are unable to initialize their connection.
    clientConnectionLost (connector, reason)
        Clients call this when their connections are lost.
class paws.core.plugins.TCPClientPlugin.TCPClientPlugin
    Bases: paws.core.plugins.plugin.PawsPlugin
    content ()
    description ()
    send_text (txt)
    start ()
    stop ()
class paws.core.plugins.TCPClientPlugin.TCPTestProtocol
    Bases: twisted.protocols.basic.LineReceiver
```



```

addCommand (cmd)
connectionLost ()
connectionMade ()
lineReceived (line)
send_lines ()

```

### paws.core.plugins.plugin module

```

class paws.core.plugins.plugin.PawsPlugin (input_names)
    Bases: object

    content ()
        PawsPlugin.content() returns a dict containing the meaningful objects contained in the plugin. The default
        implementation returns the plugin itself, keyed by 'plugin'.

    description ()
        PawsPlugin.description() returns a string documenting the functionality of the PawsPlugin. The default
        implementation returns no description.

    start ()
        PawsPlugin.start() should perform any setup required by the plugin, for instance setting up connections
        and reading files used by the plugin. The default implementation does nothing.

    stop ()
        PawsPlugin.stop() should provide a clean end for the plugin, for instance closing all connections and files
        used by the plugin. The default implementation does nothing, assumes the plugin can be cleanly terminated
        by dereferencing.

```

### paws.core.plugins.plugin\_manager module

```

class paws.core.plugins.plugin_manager.PluginManager (**kwargs)
    Bases: paws.core.treemodel.TreeSelectionModel

    Tree structure for managing pawes plugins.

    add_plugin (pgin_tag, pgin)
        Add a Plugin to the tree as a new top-level TreeItem.

    build_dict (x)
        Overloaded build_dict to handle Plugins

    get_plugin_byname (pgin_name)

    headerData (section, orientation, data_role)

    load_from_dict (pgin_dict)
        Load plugins from a dict that specifies their setup parameters.

    remove_plugin (rm_idx)
        Remove a Plugin from the tree

    staticMetaObject = <PySide.QtCore.QMetaObject object>

    write_log (msg)

```

## Module contents

`paws.core.plugins.load_plugins(path_, pkg)`

## paws.core.workflow package

### Submodules

#### paws.core.workflow.wf\_manager module

**class** `paws.core.workflow.wf_manager.WfManager(qapp_reference, plugin_manager, **kwargs)`

Bases: `paws.core.treemodel.TreeSelectionModel`

Tree structure for managing a Workflow built from paws Operations.

**add\_op** (*uri, new\_op*)

Add an Operation to the tree as a new top-level TreeItem.

**batch\_op\_stack** (*b\_itm, valid\_wf\_inputs*)

Use `b_itm.data.batch_ops()` and a list of `valid_wf_inputs` (a list of uris that are good inputs from the workflow) to build a stack (list) of lists of operations such that all Operations in the stack have their dependencies satisfied by the Operations above them.

**build\_dict** (*x*)

Overloaded `build_dict` to handle Operations. Base class method builds dicts from other data types.

**check\_wf** ()

Check the dependencies of the workflow. Ensure that all loaded operations have inputs that make sense. Return a status code and message for each of the Operations.

**columnCount** (*parent*)

Let `WfManager` have two columns, one for item tag, one for item type

**execution\_stack** ()

Build a stack (list) of lists of TreeItems, such that each TreeItem list contains a set of Operations whose dependencies are satisfied by the operations above them. For Batch or Realtime operations, the layer should be of the form `[batch_item, batch_stack]`, where `batch_item.data` is the batch controller Operation, and `batch_stack` is built from `self.batch_op_stack()`.

**finish\_thread** (*th\_idx*)

**get\_valid\_wf\_inputs** (*itm*)

Return all of the TreeModel uris of `itm` and its children which can be used as downstream inputs in the workflow.

**headerData** (*section, orientation, data\_role*)

**inputs\_tag** = 'inputs'

**is\_op\_ready** (*itm, valid\_wf\_inputs, batch\_routes=[]*)

**is\_running** ()

**load\_from\_dict** (*opman, opdict*)

Load things in to the Workflow from an `OpManager` using a dict that specifies operation setup.

**load\_inputs** (*op*)

Loads input data for an Operation from that Operation's `input_locators`. It is expected that `op.input_locator[name]` will refer to an `InputLocator`.

**locate\_input** (*il*)  
Return the data pointed to by a given InputLocator object.

**loopwait** (*interval*)  
Create an event loop to delay some time without busywaiting. Time interval is specified in milliseconds.

**next\_available\_thread** ()

**outputs\_tag** = 'outputs'

**remove\_op** (*rm\_idx*)  
Remove an Operation from the workflow tree

**run\_wf** ()

**run\_wf\_batch** (*b\_itm, stk*)  
Executes the items in the stack *stk* under the control of one Batch controller Operation

**run\_wf\_realtime** (*rt\_itm, stk*)  
Executes the workflow under the control of one Realtime controller Operation, where the realtime controller Operation is found at *rt\_itm.data*.

**run\_wf\_serial** (*stk, thd\_idx=None*)  
Serially execute the operations contained in the stack *stk*.

**set\_op\_input\_at\_uri** (*uri, val*)  
Set an op input, indicated by *uri*, to provided value. *uri* must be of the form *op\_name.inputs.input\_name*. Currently shallower uris (e.g. *op\_name.inputs*) and deeper uris (e.g. *op\_name.inputs.input\_list.0*) are not supported.

**staticMetaObject** = <PySide.QtCore.QMetaObject object>

**stop\_wf** ()

**updateOperation** (*tag, op*)

**update\_io\_deps** ()  
Remove any broken dependencies in the workflow. Should only be called after all current data have been stored in the tree.

**update\_op** (*uri, new\_op*)  
Update Operation in treeitem indicated by *uri*. It is expected that *new\_op* is a reference to the Operation stored at *uri*.

**static update\_uri\_dict** (*d, d\_new*)

**uri\_to\_dict** (*uri, data*)

**wait\_for\_thread** (*th\_idx*)  
Wait for the thread at *self.\_wf\_threads[th\_idx]* to be finished

**wait\_for\_threads** ()  
Wait for all workflow execution threads to finish

**wfdone** = <PySide.QtCore.Signal object>

**write\_log** (*msg*)

## paws.core.workflow.wf\_worker module

**class** `paws.core.workflow.wf_worker.WfWorker` (*to\_run=None, parent=None*)  
Bases: `PySide.QtCore.QObject`

Container for storing and executing parts of a workflow, to be pushed onto QtCore.QThread(s) as needed.

```
finished = <PySide.QtCore.Signal object>
opDone = <PySide.QtCore.Signal object>
staticMetaObject = <PySide.QtCore.QMetaObject object>
work ()
```

## Module contents

### Submodules

#### paws.core.listmodel module

```
class paws.core.listmodel.ListModel (input_list=[], parent=None)
    Bases: PySide.QtCore.QAbstractListModel

    Class for list management with a QAbstractListModel. Implements required virtual methods rowCount() and
    data(). Resizable ListModels must implement insertRows(), removeRows(). If a nicely labeled header is
    desired, implement headerData().

    append_item (thing)

    columnCount (parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >)

    data (idx, data_role)

    flags (idx)

    get_item (idx)

    headerData (section, orientation, data_role)

    insertRows (row, count)

    list_data ()

    removeRows (row, count, parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >)

    remove_item (row)

    rowCount (parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >)

    set_disabled (row)

    set_enabled (row)

    staticMetaObject = <PySide.QtCore.QMetaObject object>
```

#### paws.core.pawstools module

```
class paws.core.pawstools.FileSystemIterator (dirpath, regex)
    Bases: _abcoll.Iterator

    next ()

exception paws.core.pawstools.LazyCodeError (msg)
    Bases: exceptions.Exception
```

```
paws.core.pawstools.dtstr()
```

Return date and time as a string

```
paws.core.pawstools.timestr()
```

Return time as a string

## paws.core.treeitem module

```
class paws.core.treeitem.TreeItem(row, column, parent)
```

Bases: object

Container for packing objects into a TreeModel.

This is a container to facilitate data storage in a TreeModel. A TreeItem keeps references to a parent QModelIndex, and to its row and column within the QAbstractItemModel structure. The objective content of the TreeItem is stored at TreeItem.data. Every TreeItem must have a tag() for display and uri creation.

```
children_checked()
```

```
insert_child(new_child, row)
```

```
is_checked()
```

```
n_children()
```

```
remove_child(row)
```

```
set_checked(val)
```

```
set_tag(tag_in)
```

```
tag()
```

## paws.core.treemodel module

```
class paws.core.treemodel.TreeModel
```

Bases: PySide.QtCore.QAbstractItemModel

Class for tree management with a QAbstractItemModel. Implements required virtual methods index(), parent(), rowCount(). Other required virtual methods are columnCount() and data(): these should be implemented by subclassing of TreeModel. Resizable TreeModels must implement: insertRows(), removeRows(), insertColumns(), removeColumns(). If nicely labeled headers are desired, one should implement headerData().

```
auto_tag(prefix)
```

Generate the next unique tag from prefix by appending ‘\_x’ to it, where x is a minimal nonnegative integer.

```
build_dict(x)
```

Build a dict from structured data object x

```
build_uri(idx)
```

Build a URI for the TreeItem at idx by prepending its parent tags with ‘.’ as a delimiter.

```
columnCount(parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >)
```

Let TreeModels by default have one column, to display the local TreeItem’s tag.

```
data(itm_idx, data_role)
```

```
get_from_uri(uri)
```

Get from this tree the item at the given uri.

**get\_item** (*idx*)  
Just a prettier face in front of `idx.internalPointer()`

**headerData** (*section, orientation, data\_role*)

**index** (*row, col, parent*)  
Returns `QModelIndex` address of int row, int col, under `QModelIndex` parent. If a row, column, parent combination points to an invalid index, returns invalid `QModelIndex()`.

**insertRows** (*row, count, parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

**is\_good\_tag** (*testtag, parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)  
Checks for usable tags, returns a (bool,string) tuple where the bool indicates whether the tag is good, and the string provides explanation if the tag is not good.

**is\_good\_uri** (*uri*)  
Returns whether or not input uri points to an item in this tree.

**iter\_indexes** (*parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)  
Provide a list of the `QModelIndexes` held in the tree

**list\_tags** (*parent*)  
Get a list of tags for `TreeItems` under parent.

**parent** (*idx*)  
Returns `QModelIndex` of parent of item at `QModelIndex` index

**print\_tree** (*rowprefix='', parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

**removeRows** (*row, count, parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

**rowCount** (*parent=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)  
Either give the number of top-level items, or count the children of parent

**setData** (*idx, value, data\_role*)

**staticMetaObject** = `<PySide.QtCore.QMetaObject object>`

**static\_tag\_error** (*tag, err\_msg*)  
Provide a human-readable error message for bad tags.

**tree\_dataChanged** (*idx*)

**tree\_update** (*idx, x\_new*)  
Call this function to store `x_new` in the `TreeItem` at `idx` and then build/update/prune the subtree rooted at that item. Take measures to change as little as possible of the tree, since this can be a big operation and is called frequently.

**class** `paws.core.treemodel.TreeSelectionModel`  
Bases: `paws.core.treemodel.TreeModel`

**columnCount** (*parent*)  
Let `TreeSelectionModel` have two columns: one for the `TreeItem` tag, one for selection status.

**data** (*itm\_idx, data\_role*)

**flags** (*idx*)

**get\_all\_selected** (*idx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

**headerData** (*section, orientation, data\_role*)

**setData** (*idx, value, data\_role*)

**set\_all\_unselected** (*idx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

## Module contents

## Subpackages

## Submodules

```
class paws.ui.widgets.op_widget.OpWidget (op)
    Bases: PySide.QtGui.QWidget

    paintEvent (evnt)

    staticMetaObject = <PySide.QtCore.QMetaObject object>
```

```
class paws.ui.widgets.pif_widget.PifWidget (itm)
    Bases: PySide.QtGui.QTextEdit

    print_comp (itm, indent)
    print_id (id_, indent)
    print_matrix (itm, indent)
    print_pif (itm, indent)
    print_procstep (itm, indent)
    print_prop (itm, indent)
    print_qty (itm, indent)
    print_scalar (itm, indent)
    print_src (itm, indent)
    print_value (itm, indent)
    print_vector (itm, indent)

    staticMetaObject = <PySide.QtCore.QMetaObject object>
```

```
paws.ui.widgets.text_widgets.display_text(itm, indent='&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;')  
paws.ui.widgets.text_widgets.display_text_fast(itm,                               in-  
                                                dent='&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;&nbsp;')
```

## Module contents

widgets package

## Submodules

### paws.ui.data\_viewer module

`paws.ui.data_viewer.display_item(itm, uri, qlayout, logmethod=None)`

### paws.ui.input\_loader module

`class paws.ui.input_loader.InputLoader(name, src, trmod=None, parent=None)`

Bases: `object`

This class controls the `input_loader.ui` to select inputs from various sources.

`add_items()`

`add_value(val)`

`lock_list_toggle()`

`rm_selection()`

`set_list_toggle()`

`setup_ui()`

`unlock_list_toggle()`

`unset_list_toggle()`

### paws.ui.opuiman module

`class paws.ui.opuiman.OpUiManager(opman)`

Bases: `object`

`build_source()`

`clear_op()`

`edit_op()`

`enable_ops()`

`get_op(idx)`

`set_op(op)`

Set up ui elements around existing input op

`setup_ui()`



### paws.ui.plotmaker\_mpl module

```
paws.ui.plotmaker_mpl.array_plot_1d(data_in)
paws.ui.plotmaker_mpl.array_plot_2d(data_in)
paws.ui.plotmaker_mpl.mpl_array_plot_1d(data_in)
paws.ui.plotmaker_mpl.mpl_array_plot_2d(data_in)
paws.ui.plotmaker_mpl.plot_mpl_fig(fig_in)
```

### paws.ui.plotmaker\_pqg module

```
paws.ui.plotmaker_pqg.array_plot_1d(data_in)
paws.ui.plotmaker_pqg.array_plot_2d(data_in)
paws.ui.plotmaker_pqg.plot_mpl_fig(fig_in)
paws.ui.plotmaker_pqg.pqg_array_plot_1d(data_in)
paws.ui.plotmaker_pqg.pqg_array_plot_2d(data_in)
```

### paws.ui.pluguiman module

```
class paws.ui.pluguiman.PluginListModel(input_list=[], parent=None)
    Bases: paws.core.listmodel.ListModel
    Just a ListModel with overloaded headerData
    headerData (section, orientation, data_role)
    staticMetaObject = <PySide.QtCore.QMetaObject object>
class paws.ui.pluguiman.PluginUiManager(plugman)
    Bases: object
    build_input()
    clear_input()
    fetch_data(name)
    fetch_from_input_ui(ui)
    load_plugin()
        Package the finished Plugin, ship to self.plugman
    reset_input_headers()
    reset_type_widget(name, row, src=None)
    reset_val_widget(name, row, src=None, tp=None)
    set_input(name, src_ui=None)
    set_plugin(pgin)
    setup_ui()
    start_plugin(idx)
```

**stop\_plugin()**  
remove the selected plugin from self.pluginman

## paws.ui.ui\_manager module

**class** `paws.ui.ui_manager.UiManager` (*opman, wfman, plugman*)  
Bases: `object`

Stores a reference to a QMainWindow, performs operations on it

**connect\_actions()**  
Set up the works for buttons and menu items

**display\_item** (*idx*)  
Display selected item from the workflow tree in image\_viewer

**edit\_ops** (*item\_indx=None*)  
interact with user to enable existing Operations and edit or develop new Operations

**edit\_wf** (*trmod, itm\_idx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)  
Interact with user to edit the workflow. Pass in a TreeModel and index to open the editor with the item at that index loaded.

**final\_setup()**

**make\_title()**  
Display the paws logo in the image viewer

**msg\_board\_log** (*msg, timestamp=<function timestr>*)  
Print timestamped message to msg board

**start\_load\_ui()**

**start\_plugins\_ui()**

**start\_save\_ui()**

**start\_wf\_editor** (*trmod=None, indx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)  
Create a WfUiManager (QMainWindow), return it

**toggle\_run\_wf()**

## paws.ui.uitools module

Configuration flags, widgets, and functions for paws gui control.

`paws.ui.uitools.bigtext_widget` (*text*)

`paws.ui.uitools.hdr_widget` (*text*)

`paws.ui.uitools.load_path` (*ui, idx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

`paws.ui.uitools.message_ui` (*parent=None*)

`paws.ui.uitools.name_widget` (*name*)

`paws.ui.uitools.r_hdr_widget` (*text*)

`paws.ui.uitools.save_path` (*ui, idx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >, oldidx=<PySide.QtCore.QModelIndex(-1, -1, 0x0, QObject(0x0)) >*)

`paws.ui.uitools.save_to_file` (*d, filename*)  
Save the items in dict d as YAML in filename

```

paws.ui.uitools.smalltext_widget (text)
paws.ui.uitools.src_selection_widget ()
paws.ui.uitools.start_load_ui (uiman)
    Start a modal window dialog to load a previously saved workflow
paws.ui.uitools.start_save_ui (uiman)
    Start a modal window dialog to choose a save destination for the current workflow
paws.ui.uitools.stop_load_ui (ui, uiman)
paws.ui.uitools.stop_save_ui (ui, uiman)
paws.ui.uitools.text_widget (text)
    Produce a Read-only Center-aligned QtGui.QLineEdit from input text.
paws.ui.uitools.toggle_expand (trview, idx)
paws.ui.uitools.toggle_load_button (ui, txt)
paws.ui.uitools.toggle_save_button (ui, txt)
paws.ui.uitools.type_selection_widget (src, widg=None)

```

## paws.ui.wfuiman module

```

class paws.ui.wfuiman.WfUiManager (wfman, opman, plugman)
    Bases: object

    build_input_locator (name, ui=None)
        Create an InputLocator for named input. If a ui is provided, it should be an input_loader.ui, and this should
        load data from that ui.

    build_io ()

    clear_io ()

    create_op (op)
        Instantiate op, call self.set_op()

    fetch_data (name)

    fetch_from_input_ui (ui)

    get_op (trmod, itm_idx)

    load_op ()
        Package the finished Operation, ship to self.wfman

    reset_input_headers ()

    reset_output_headers ()

    reset_type_widget (name, row, src=None)

    reset_val_widget (name, row, src=None, tp=None)

    rm_op ()
        remove the selected operation from the workflow

    set_input (name, src_ui=None)
        Call build_input_locator to package an InputLocator for named input. Store it in
        self.op.input_locator[name]. This is called on all inputs when an Operation is loaded, so it should
        not alter an identical already-loaded input.

```

**set\_op** (*op, uri*)

Set up ui elements around existing input op

**setup\_ui** ()

**Module contents**

**Module contents**

## INDICES AND TABLES

- genindex
- modindex
- search



## PYTHON MODULE INDEX

### p

paws, 40  
 paws.core, 35  
 paws.core.listmodel, 32  
 paws.core.operations, 27  
 paws.core.operations.COMM, 9  
 paws.core.operations.COMM.TestTCP, 9  
 paws.core.operations.DISPLAY, 10  
 paws.core.operations.DISPLAY.simple\_plots, 10  
 paws.core.operations.DISPLAY.zip, 10  
 paws.core.operations.EXECUTION, 11  
 paws.core.operations.EXECUTION.batch\_execution, 10  
 paws.core.operations.EXECUTION.realtime\_execution, 11  
 paws.core.operations.INPUT, 13  
 paws.core.operations.INPUT.fabio\_input, 12  
 paws.core.operations.INPUT.pillow\_tif\_input, 12  
 paws.core.operations.INPUT.read\_csv, 13  
 paws.core.operations.INPUT.WXDIFF, 12  
 paws.core.operations.INPUT.WXDIFF.read\_wxdiff\_files, 12  
 paws.core.operations.op\_manager, 25  
 paws.core.operations.operation, 26  
 paws.core.operations.optools, 27  
 paws.core.operations.OUTPUT, 14  
 paws.core.operations.OUTPUT.PIF, 14  
 paws.core.operations.OUTPUT.PIF.citrination\_shipment, 13  
 paws.core.operations.OUTPUT.write\_csv, 14  
 paws.core.operations.PACKAGING, 15  
 paws.core.operations.PACKAGING.data\_packing, 15  
 paws.core.operations.PACKAGING.PIF, 15  
 paws.core.operations.PACKAGING.PIF.pif\_saxs, 14  
 paws.core.operations.PACKAGING.PIF.pif\_ternary, 15  
 paws.core.operations.PROCESSING, 25  
 paws.core.operations.PROCESSING.arithmetic, 21  
 paws.core.operations.PROCESSING.calib\_and\_reduce, 22  
 paws.core.operations.PROCESSING.HiTp\_feature\_extraction, 20  
 paws.core.operations.PROCESSING.mirror, 22  
 paws.core.operations.PROCESSING.rotation, 23  
 paws.core.operations.PROCESSING.SAXS, 20  
 paws.core.operations.PROCESSING.SAXS.background\_subtraction, 16  
 paws.core.operations.PROCESSING.SAXS.saxs, 16  
 paws.core.operations.PROCESSING.SAXS.saxs\_features, 20  
 paws.core.operations.PROCESSING.smoothing, 23  
 paws.core.operations.PROCESSING.VoigtPeakFit, 21  
 paws.core.operations.TESTS, 25  
 paws.core.operations.TESTS.identity, 25  
 paws.core.pawstools, 32  
 paws.core.plugins, 30  
 paws.core.plugins.CitrinationPlugin, 27  
 paws.core.plugins.plugin, 29  
 paws.core.plugins.plugin\_manager, 29  
 paws.core.plugins.SpecClientPlugin, 28  
 paws.core.plugins.TCPClientPlugin, 28  
 paws.core.treeitem, 33  
 paws.core.treemodel, 33  
 paws.core.workflow, 32  
 paws.core.workflow.wf\_manager, 30  
 paws.core.workflow.wf\_worker, 31  
 paws.ui, 40  
 paws.ui.data\_viewer, 36  
 paws.ui.input\_loader, 36  
 paws.ui.opuiman, 36  
 paws.ui.plotmaker\_mpl, 37

`paws.ui.plotmaker_pgg`, [37](#)  
`paws.ui.pluguiman`, [37](#)  
`paws.ui.ui_manager`, [38](#)  
`paws.ui.uitools`, [38](#)  
`paws.ui.wfuiman`, [39](#)  
`paws.ui.widgets`, [36](#)  
`paws.ui.widgets.op_widget`, [35](#)  
`paws.ui.widgets.pif_widget`, [35](#)  
`paws.ui.widgets.text_widgets`, [35](#)