
paws Documentation

Release 0.4.0

Lenson A. Pellouchoud

Feb 21, 2017

CONTENTS

1	Introduction	3
2	Quick Start	5
3	Installation	7
3.1	Installing with pip	7
3.2	Downloading Source	7
3.3	Testing	7
4	API Documentation	9
4.1	paws package	9
5	Indices and tables	11
	Python Module Index	13

Contents:

INTRODUCTION

The `paws` package aims to provide a fast and lean platform for building and executing workflows for data processing. It was originally developed to perform analysis of diffraction images for research purposes at SLAC/SSRL. At the core of `paws` is a workflow engine that uses a library of operations to crunch through data and expose select results while attempting to minimize resource consumption.

`paws` is currently written in Python, based on Qt via the PySide bindings. Internally, `paws` keeps track of data in Qt-based tree models, which can be controlled either directly (through the `paws` api) or through a gui (employing the Qt model-view framework).

`paws` also provides an interface to `xi-cam`, a synchrotron x-ray diffraction data analysis package written by the CAMERA Institute and Pandolfi, et al at the Lawrence Berkeley National Lab.

Some the core goals of `paws`:

- Eliminate redundant development efforts
- Streamline and standardize routine data analysis
- Simplify data storage and provide large-scale analysis
- Perform data analysis in real time for results-driven feedback

The `paws` developers would love to hear from you if you have wisdom, thoughts, haikus, bugs, artwork, or suggestions. Limericks are also welcome. Get in touch with us at `paws-developers@slac.stanford.edu`.

QUICK START

Minimal and usually-effective installation instructions.

Here is a reference to the *[brief introduction](#)*.

This chapter is for setting up `paws` quickly in an environment that is prepared to install Python packages with `pip`.

INSTALLATION

Here are instructions for installing `paws` from PyPI, or downloading and testing the `paws` source code.

Installing with pip

Instructions will go here for installing `paws` using the Python package installer `pip` (currently not implemented).

Downloading Source

The source code for `paws` is hosted on `github`. Clone the repository from `https://github.com/slaclab/paws.git`. You should then be able to run `paws` by invoking `python main.py` from the root directory.

Testing

`paws` comes with a tests that can be used to ensure the platform runs as expected. After *downloading the source*, invoke `python -m unittest discover` from the root directory.

API DOCUMENTATION

This is the complete auto-generated documentation of the `paws` package, made with sphinx-apidoc.

paws package

Subpackages

paws.api package

Module contents

Module defining the API for `paws`

`paws.api.core_app(app_args=[])`

Return a reference to a new `QCoreApplication` or a currently running `QApplication`.

Input arguments are passed to the `QApplication` constructor. If a `RuntimeError` is thrown, it is assumed that a `QApplication` is already running, and an attempt is made to return a reference to that `QApplication`. If that fails, this returns `None`.

Parameters `app_args` – arguments to pass to the `QApplication` constructor

Returns reference to a new or existing `QCoreApplication`

Return type `PySide.QtCore.QCoreApplication` or `None`

`paws.api.start(app_args=[])`

Instantiate an `Operation Manager`, a `Workflow Manager`, and a `Plugin Manager`. Return references to them.

`paws.api.start()` calls `paws.api.core_app()`, then sets up and returns references to a `paws Workflow Manager` (`paws.api.workflow_manager`), `Operation Manager` (`paws.api.op_manager`), and `Plugin Manager` (`paws.api.plugin_manager`).

Parameters `app_args` – arguments to pass to the `QApplication` constructor

Returns references to `paws operation manager`, `workflow manager`, and `plugin manager`

Return type tuple of `paws.core.operations.op_manager.OpManager`,
`paws.core.workflow.wf_manager.WfManager`, `paws.core.plugins.plugin_manager.PluginManager`

Module contents

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

p

paws, [9](#)

paws.api, [9](#)